



Determination of the Formula of a Hydrate

Procedure:-


- ① use tongs all time:-




- ② Heat up the crucible:- then it need to cool



- ③ Get the mass of crucible and lid using Analytical Balance:-



- ④ weight 0.5 approx using weight paper! hydrate ~~etc~~




- ⑤ Put hydrate to crucible:- then get mass of both.



- ⑥ Heat the crucible for 10-15 minutes or to maintain temp
let steam escape




- ⑦ Transfer crucible to desiccator with lid.



let it cool.


- ⑧ record the mass then again reheat and cool record the mass. Repeat it till final mass



cool



- ⑨ Add drops of water.



notice the change

10: Calculate no. of moles of anhydrous and water moles from your hydrate sample.