



Determination of the molecular weight of An Unknown acid.


- (1) Prepare a standard solution of $\text{H}_2\text{C}_2\text{O}_4 \cdot 2\text{H}_2\text{O}$
(2) weight approx 1.5g of acid $\pm 0.02\text{g}$, and weight dry beaker




- (3) Add oxalic acid to beaker, and weight it.




- (4) add crystal to distil water



- (5) Transfer solution to 250 ml flask




- (6) Wash beaker and stirring rod with distilled water




⑦

mix gently, and add water with dropper




⑧

transfer oxalic acid 25 ml to flask




⑤

add 3 drops of phenolphthalein indicator



⑥

get reading of initial and final.



repeat the cycle



until result. then do

⑦


weight the unknown acid calculations



± 0.002

8.

Transfer sample into 250 ml Erlenmeyer flask



9.7

Add drops of phenolphthalein



10.7 Do calculations!